

CNSeq

Advanced Chromosomal Analysis with
Low-pass whole genome sequencing





CNSeq: Low-pass Whole Genome Sequencing

CNSeq revolutionizes chromosomal analysis. Leveraging NGS technology, it detects Copy Number Variants (CNVs) with unmatched accuracy compared to traditional cytogenetic methods.

The key to CNSeq's power lies in its strategic use of low-pass whole genome sequencing, reading each base a few times. This focuses on identifying large-scale CNVs uniformly across the genome for faster analysis and potentially deeper insights into chromosomal abnormalities.

It empowers clinicians with the ability to explore insights in detail, enabling a comprehensive grasp of structural variations in the genome with unbiased CNV calling across the genome.



CNSeq



CNSeqHD

Key Features



Increased CNV callers for improved detection reliability



Validated for both prenatal & pediatric settings



Advanced pipeline adjusts for maternal cell contamination in calculating mosaicism and aneuploidies



Turnaround Time
21 Days




Sample Requirement
Amniotic Fluid, CVS, POC, PVB & gDNA



Going beyond traditional array with CNSeq

Technique	Karyotype	Prenatal CMA	750K	CNSeq	CNSeqHD
Coverage	>5Mb	Probe-based	Probe-based	Complete genome coverage	Complete genome coverage
Resolution	>5Mb	>500Kb	>100Kb	500Kb	>50Kb
DNA quantity requirement	Cell Based	250 ng	250 ng	50 ng	50 ng
Low-level Mosaicism sensitivity	10%	>20%	>20%	>20%	>20%
Abnormalities detected	Aneuploidies, larger CNVs, Structural Abnormalities, Balanced Translocation	Aneuploidies, CNVs, UPDs, Polyploidy, AOH/LOH	Aneuploidies, CNVs, UPDs, Polyploidy, AOH	All Aneuploidies, CNVs & Structural Variants	All Aneuploidies, CNVs & Structural Variants






CNSeq & CNSeqHD: Validation Data

Coriell Controls	Detected Events – CNVs & Aneuploidies (bp)							Whole Chr
	10 – 100K	100 – 250K	250 – 500K	500K – 1M	1M – 5M	5M – 30M	>30M	
CNSeq								
NA50324								
NA21699								
NA04375								
NA04409								
NA03623								
NA11385								
CNSeqHD								
NA21886								
NA04409								
NA12878								
PCE-2								
PCE-3								
NA21699								

● Gain; ■ Mosaic Gain; ▲ Loss; ◆ Polymorphic event undetected; ■ Borderline gain event undetected

Unique CNVs picked up by CNSeq

Wolf-Hirschhorn Syndrome - Pathogenic loss



Unique CNVs picked up by CNSeqHD

Congenital Heart Defect – 10kb Deletion Pathogenic Copy Number Loss



Proprietary Platforms

strandomics

Proprietary platform for unparalleled Variant Interpretation & Reporting

strandngs

State-of-the-art platform for Best in class Sequencing

illumina
NovaSeq X Plus

World's highest throughput Sequencer for fastest report delivery

Additional Testing Options

- ExomePlus: Exome Sequencing
- MaatriSeq: NIPS
- AneuXpress: Quantitative Fluorescent PCR
- Sanger Based Testing
- Couple Carrier Screening
- Whole Genome Sequencing
- Preimplantation Genetic Testing: PGT – Aneuploidy

24⁺

YEARS OF
EXPERIENCE

80,000+
Genetic Tests
Reported

500+ Projects
Executed for
Genomics
Majors Globally


Presence in
20+ Countries


SCAN TO KNOW MORE




CHANNEL
PARTNER

strand 

 Reference Sdn Bhd (1278003-P) Wisma MEDINICS, No.2, Jalan Astaka U8/88B, Seksyen U8, Bukit Jelutong, 40150 Shah Alam, Selangor, Malaysia.

 +6019 400 3352

 support@refmedlab.com

 www.refmedlab.com